




RIGID Head and Neck palm model Ref: 101-3+tank

Arterio-venous palm model in rigid silicone with tank and pump

Arterio - Venous model



Ref 101-3+tank = Delivery of a complete set ready to work in a closed water circuit.

A horizontal tank with an aquarium continuous flow pump and the prolonged rigid silicone model, to work in a closed water circuit.

Approx external dimensions of rigid silicone bloc: 16.5cm x 13cm x 1.60 cm

This model simulates part of an arterio - venous malformation and is particularly indicated for embolization trainings.