



EXAMPLE OF A VASCULAR BODY MODEL


Models shown here under: Thorax T-S-N-002+, Left and Right ARMS, Abdominal A-S-N-002 + Left and Right Femorals and Tibials.

This version is shown with the tube inside the boxes. No water can be inserted into the boxes with this version!

Should the boxes be filled to mimic as close as possible the human body, the tubes need to be outside the boxes and Elastrat informed.

The client's version can be chosen within our large range of head, thoracic and abdominal models = shown on our website

www.elastrat.com. A continuous flow aquarium pump works in a closed water circuit. A pulsatile pump can also be plugged.



ELASTRAT's models received the **SWISS MADE LABEL**, official SWISS Certification for many years.