




RIGID Head and Neck palm model Ref: 101-1+

Arterio-venous palm model available in a transparent rigid silicon



This model simulates part of an arterio-venous malformation and is particularly indicated for the training of embolization products, as the Onyx would be.


Dimensions of the rigid silicon palm model: 15 x 11 x 2 cm.

By using similar products as the Onyx, each model can be cleaned four times (4x) by Elastrat.



RIGID Head and Neck palm model Ref: 101-1+

Version with a 3 liter **plexiglas** tank and an aquarium pump.
Choice between a 3 lt Plexiglas tank here under or a plastic tank of 2 lt.



Execution with the 2 liter **plastic** tank:

