





Screw Introducers of 6 French, 12F, 18F and 24F:

Elastrat's Introducer let you introduce your devices up to 24F = 8mm




Example of an Introducer used with an abdominal model Ref ASN-002+ (right side of picture).



This screw introducer goes from wide open to completely closed with its soft silicone integrated ring.
 Piece 1 = stopper to be introduced uniquely if you want the device being completely closed.
 Piece 2 = the exchangeable white silicon piece. Available in 6 French = 2mm, 12F = 4mm, 18F = 6mm, 24F = 8mm.

Please do release (unscrew) the introducer when not in use. The inside soft silicone seal has to come back to no pressure otherwise the white silicone seal will be crushed with time.


Picture to show you the white seal (2) to exchange it if needed.





Example of Screw Introducers with Internal Diameter of: 6 French, 12F, 18F and 24F:

Different stoppers (1) and their respective soft silicone pieces (2).
The stopper goes perfectly into the soft silicone piece to obtain the seal.



Have a good use with the models and their screw introducers.

The Elastrat Team.